


PHOENIX ANESTHESIA


TIVA Pump




induction screen



maintenance screen



bolus mode selection screen



history screen



This device complies with MDD 93/42/CEE

FOURES SAS, Infusion Department

Voie romaine, parc Remora, 33173 GRADIGNAN FRANCE

Tél. (33) 05 56 89 26 60 - Fax (33) 05 56 89 54 42 - sasfoures@foures.fr

Specifications

20 drug protocols including

- drug name, drug concentration
- induction dose, induction time
- maintenance rate
- bolus dose, bolus time

5 infusion scheme protocols including

- infused drug (among the 20 selectable drugs)
- maintenance rate and time for each stage (3 max)

induction/bolus rate	: up to 1200 ml/h
drug concentration	: 1µg/ml - 99.9 mg/ml
patient weight	: 0.1 - 200.0 kg
occlusion limits	: 225/450/675/900 mmHg (0.3/0.6/0.9/1.2 bar)

Master-Slave quick programmation mode
RS 232 serial port for printer
History printing without the need of a computer

Power supply : 230 V - 50 Hz
internal battery - 6 hrs autonomy

Dimensions : 272 x 130 x 170 mm
Weight : 3.3 Kg
ABS case
Options : thermal printer, Master-Slave cable

Le fabricant se réserve le droit de modifier à tout moment la gamme PHOENIX sans avis préalable. Cette plaquette ne pourra être considérée comme un document contractuel.

FOURÈS S.A.

INFUSION DEPARTMENT




Total Intravenous Anesthesia is now as **easy** as you can imagine

PHOENIX ANESTHESIA

TIVA Syringe Pump

Ref : 905035529/A

Total Intravenous Anesthesia is now as easy as you can imagine

- Fast induction and bolus injections
 - Ergonomic and strong design
 - Rate changes without stopping infusion
 - Clear and intuitive front panel
 - Fast and easy starting
 - Up to 20 drug protocols including induction and maintenance parameters for fast and easy infusion management
 - Up to 5 scheme protocols automatizing drug administration according to doctor's infusion profiles
 - Hands-free or hands-on bolus
 - High performances syringe pump based on up-to-date technology
 - Full events history (up to 1600 events)
 - Back pressure permanent monitoring
 - RS 232 serial port for direct printing of pump history
- 



FOURÈS S.A.
INFUSION DEPARTMENT